

# Red Rock Medical Supply

Toll Free: 888-757-3343

E-Mail: [RedRockMedical@hotmail.com](mailto:RedRockMedical@hotmail.com)

## Eclipse Oxygen System Accessories

# 2801-7082

\$349



### Power Cartridge (battery)

- ◆ 195 Watt-hour rechargeable battery
- ◆ Extend excursion time by up to 4.4 hours (exact excursion time will vary according to oxygen flow setting)
- ◆ No tools needed to slide in or out of the Eclipse
- ◆ Automatically recharges when Eclipse is connected to either AC or DC power

# 2801-7112

\$150



### Desk Top Charger

- ◆ Great item to have for maximum versatility and mobility wherever you go
- ◆ Recharges Power Cartridge (battery) in less than 3 hours
- ◆ Can be used anywhere in the world with 100 VAC to 240 VAC, and from 50-60 Hz electrical power.
- ◆ Comes with detachable North American Power Cord
- ◆ Alternate power cords for other countries sold separately

# 2801-7086

\$145



### AC Power Supply

- ◆ Stop lugging your power supply from room to room
- ◆ Keep an extra power supply upstairs or next to your bed
- ◆ Plugs into a standard household outlet anywhere in the world with 100 VAC to 240 VAC, with 50 Hz to 60 Hz
- ◆ Comes with detachable North American Power Cord
- ◆ Alternate power cords for other countries sold separately

# 2801-7088

\$136



### DC Power Supply

- ◆ Eliminate the need to move your DC power supply from car to car
- ◆ Keep a DC power supply in each of your vehicles
- ◆ Designed to plug into any vehicle accessory outlet including auto, truck, and boat

**To Order call toll free 888-757-3343**

Payment by check, money order, Visa, MC, Discover, Amex