

# Red Rock Medical Supply

Toll Free: 888-757-3343

Email: [RedRockMedical@hotmail.com](mailto:RedRockMedical@hotmail.com)

## EverFlo™ Stationary Oxygen Concentrator

by Respironics

The quiet, lightweight oxygen solution

### Features

- ◆ **Lightweight** and easy to transport at 31 pounds
- ◆ **Compact**, ergonomic design takes up less space, looks less like a medical device
- ◆ **Quiet**, you will need to check the power light to be sure it is running
- ◆ **Uses less power** and produces less heat
- ◆ **Recessed flow meter** helps reduce the chance of accidental breakage



### Specifications

Input Voltage	120 VAC (+/- 10%)
Input Frequency	60 Hz
Average Power Consumption	350 Watts
Oxygen Concentration (at 5 LPM)	93% (+/- 3%)
Liter Flow	0.5 - 5.0 liters per minute
Weight	31 lbs
Sound Level ( <b>Everflo Standard</b> )	45 dba typical
Sound Level ( <b>Everflo Q</b> )	40 dba typical
Dimensions (HxWxD)	23 in H x 15 in W x 9.5 in D
Operating Temperature	55 to 90 degrees F
Operating Humidity	up to 95% relative humidity
Outlet Pressure	5.5 PSI
Operating Altitude	0 to 7500 feet

### Price

#2702-1020000 EverFlo Standard Oxygen Concentrator	<b>\$ 469</b>
#2702-1020015 EverFlo Q Ultra Quiet Oxygen Concentrator	<b>\$ 595</b>

**Red Rock pays the shipping!**

**To Order call toll free 888-757-3343**

Offer good while supplies last.

Payment by check, money order, Visa, MC, Discover, Amex