

Red Rock Medical Supply

Toll Free: 888-757-3343

Web Site: RedRockMedicalSupply.com

Compressor Nebulizer Systems by Respironics

#2702-RRD481

MiniElite™

\$ 169



System complete with compressor, reusable nebulizer and tubing, AC & DC adapters, rechargeable battery, air filter pack, carry case, and instructional information

- ✦ Quiet, small, and only 1.4 lb. with battery
- ✦ Treatment times in as little as seven minutes
- ✦ Up to 90 minutes of treatment before battery needs recharge
- ✦ Easy to use, transport, and maintain
- ✦ Nebulizes most commonly prescribed solutions & suspensions
- ✦ Three-year warranty on compressor

#2702-RRD491

MicroElite™

\$ 189



System complete with compressor, reusable nebulizer handset, AC & DC adapters, rechargeable battery, four filters, carry case, and instructional DVD

- ✦ Quiet, small, and only 1 lb. with battery
- ✦ Treatment times in as little as eight minutes
- ✦ Up to 60 minutes of treatment before battery needs recharge
- ✦ Easy to use, transport, and maintain
- ✦ Nebulizes most commonly prescribed solutions & suspensions
- ✦ Three-year warranty

Specifications

	2702-RRD491	2702-RRD481
Weight with battery	1.0 lbs.	1.4 lbs.
Size with battery	3" x 2.7" x 2.5"	4.5" x 3.8" x 2.5"
Nebulization Capacity	6 ml.	8 ml.
Treatment Time	8 minutes for 3 ml.	7 minutes for 3 ml.
Operating Pressure	9.4 PSI	11.6 PSI
Operating Flow Rate	2 LPM	3.5 LPM
Battery Discharge Time	60 minutes	90 minutes

To Order call toll free **888-757-3343**

Payment by check, money order, Visa, MC, Discover, Amex