

Red Rock Medical Supply

Toll Free: 888-757-3343

Web Site: www.RedRockMedicalSupply.com

Oxygen System Accessories

1302-DHO6M9

\$69



Six Pak Carrier

- ◆ Organize and carry extra cylinders easily and safely
- ◆ Holds 6 M9, M6, M4, or M2 cylinders
- ◆ Compact 11" W x 15" L x 14.5" H
- ◆ Lightweight at only 5 pounds
- ◆ Durable chrome plated steel
- ◆ Oxygen cylinder shown is sold separately

1302-HC210

\$89



Double Cylinder Cart

- ◆ Double your excursion time with two cylinder capacity
- ◆ Safely and easily transport two E or D size cylinders
- ◆ Black non marring wheels do not show dirt
- ◆ Adjustable handle for right or left hand use
- ◆ Handle height adjustment from 32" to 42"
- ◆ Front lifting handle included
- ◆ Stable 12" x 13" footprint
- ◆ Weighs just 12 pounds
- ◆ Oxygen cylinders shown are sold separately

2803-PRO2

\$595



Concentrator Analyzer

- ◆ Oxygen and air taste and smell the same!
- ◆ Most concentrators continue to run even when they no longer concentrator air into oxygen.
- ◆ Think you are just having a bad sick day? Chances are very good that it is your concentrator that is sick.
- ◆ Attaches like a cannula to easily and accurately determine the % oxygen concentration of any oxygen device.
- ◆ Powered by a single 9 volt battery, included.

To Order call toll free **888-757-3343**

Payment by check, money order, Visa, MC, Discover, Amex